

Pineapple[®] solo MH

US USER MANUAL

PRO00096_US
Revision 1.0



Pineapple

manufactured by **medstrom⁺**

This page is intentionally blank.

Contents

1.	Notifications Used in the Manual	5
2.	Identification of Labels Affixed to Bed	5
3.	Intended Use of the Bed	8
4.	General Information	8
5.	Precautions for Use	8
5.1.	General Precautions	8
5.2.	Bed Frame Compatibility	10
5.3.	Contraindications	10
5.4.	Potential Risk of Entrapment	11
5.5.	Installation and Commissioning	12
6.	Transport and Storage	14
6.1.	Shutdown Procedure	14
7.	Bed Features	14
7.1.	Dimensions	16
7.2.	Technical Specifications	17
8.	Instructions for Use	18
8.1.	Wireless Handset	18
8.2.	Recharging the Battery	20
8.3.	Variable Bed Height Function	20
8.4.	Sleep Deck Articulation Functions	22
8.5.	Braking System	26
8.6.	Steering System - AutoSteer™	27
8.7.	Moving the Bed	28
9.	Accessories:	29
9.1.	Mental Health Fixed Castor Kit	29
10.	Operational Maintenance	29
10.1.	Cleaning Guidelines	29
10.2.	Sources of Dirt and Contamination	30
10.3.	Recommended Cleaning Procedures	30
11.	Electromagnetic Compatibility (EMC):	33
12.	Warranty, Servicing and Maintenance	36
13.	Service Checklist Form	37
14.	Manufacturer's Details	38




This User Manual contains safety information and operating instructions.
Please retain for future reference.



1. Notifications Used in the Manual

	Warning: There is a risk to safety. Please read this safety instruction and note any warning labels on the bed.
	Caution: There is a risk of damage to equipment. Please read this instruction.
Important:	Important information that provides further details for a procedure or condition.
Note:	A note that provides further explanation or advice.




2. Identification of Labels Affixed to Bed

 Improved Patient Outcomes  2 Cygnus Court, Beverley Road, Pegasus Business Park, Castle Donington, Derby, DE74 2SA Made in Great Britain T: +44 (0)345 371 1717 E: info@medstrom.co.uk www.medstrom.com   IPX4   =425.5lb  =568.8lb Electrical motor duty cycle 10% 2 min ON / 18 min OFF 100-240 V AC, 50/60 Hz 3.9A MAX. Voltage and Auxiliary output current 24V ~ 5A	
Label indicating the location of a possible trapping hazard:  Risk of trapping feet.	
Label indicating the location of a possible trapping hazard:  Risk of trapping fingers/hands.	
 X X X X X X X X  (01) 0XXXXXXXXXXXXX (11) XXXXXX (21) XXXXXXXXXXXXXXXX  XXXXXXXX  < Bed variant description >  YYMMDD	

Continued...

 Earth Reference Label	
Scan & Watch the Pineapple Training Video PRO00089, Rev1, CRN000338  Label for User Video	
Pineapple manufactured by medstrom⁺ Pineapple SOLO MH Branding Label	
 Brake   AutoSteer  Brake Pedal Label	
 AutoSteer   FreeMove  AutoSteer Pedal Label	
PASS ELECTRICAL TEST S/N: TEST DATE: NEXT TEST: TESTED BY: Electrical Safety Test Label. Demonstrates that the bed has been tested.	
   Label indicating where to recharge the bed, plus a dynamic QR code to access online instructions for the use of the bed.	

Continued...

 Recommended patient size for the bed.	
 Manual CPR label.	

Reference Label Symbols:

	Waste Electrical and Electronic Equipment (WEEE) Directive	Waste Electrical and Electronic Equipment must not be disposed of through general waste streams. Please contact an authorised disposal facility or Medstrom Ltd.
	Refer to User Manual	Refer to this manual before operating.
	Operating Instructions	Identifies the user manual that contains the safety information and operating instructions for the equipment.
	Protective Earth (Ground)	Terminal is intended for connection to an external conductor for protection against electrical shock in case of a fault, or is the terminal of a protective earth (ground) electrode.
	Equipotentiality (Equipotential Earth)	Identifies terminals which, when connected together, bring the various parts of the equipment to the same potential.
	Type B applied part	Identifies a type B applied part complying with IEC 60601-1. <i>Classification of protection against electrical shock.</i>
	Caution	There is a risk of damage to equipment.
	Ingress Protection (IP) Rating	IPX4: Water splashing against the equipment from any direction will have no harmful effect.
	Medical Device	Identifies the equipment as a medical device.

Continued...

	Maximum Patient Weight	The maximum weight of patient that the equipment is rated for (425.5 lb).
	Maximum Safe Working Load	The maximum safe working load (weight) (568.8 lb).
	Manufacturer	The manufacturer of the equipment (i.e. Medstrom) and their address.
	Reference	Manufacturer's reference catalogue number/identifier, i.e. SOLO MH.
	Date of Manufacture	The date that the equipment was manufactured.

3. Intended Use of the Bed

This bed is designed for long-duration use in psychiatric and mental health environments located in hospitals, nursing homes, rehabilitation centres and elderly care facilities. It is designed to aid comfortable sleep and provide smooth transfers into and out of bed. The bed ensures the highest level of care, for patients and staff, with a safe working height. The bed is designed for IEC 60601-2-52 Application Environments 2, 3 and 5, where medical care or medical supervision and control is provided, and where electro-medical devices are provided to help maintain and/or improve the condition of patients.

4. General Information

Important:

Read this user manual carefully before use and retain for future reference.

5. Precautions for Use

Important:

Read this user manual carefully before first use of the bed and any subsequent bed maintenance.

5.1. General Precautions

The following is a list of general precautions which should be followed when using the bed:



Warnings:

The bed is an electrical device with a risk of electric shock. Personnel using the bed must be informed and trained to potential risks related to electrical appliances.

It is essential to know how the product works to achieve the desired result during handling.

It is essential to comply with the use and recommendations described in this manual to ensure the security and integrity of the different users.

- The use and handling of a medical bed can cause injury if due precautions are not taken. It is essential that any person handling the bed is authorised, trained and made aware of the consequences of all actions and how not to cause involuntary movements.
- Incompatible mattresses can create hazards. See section “7.1. Dimensions” for dimensions and recommendations.
- If using the bed with other medical equipment not supplied by Medstrom, e.g. mattresses, the compatibility of the equipment must be assessed for potential risks.
- Do not use the bed for patients not meeting the patient recommendations given in section “7.1. Dimensions”.
- Do not allow sitting on the bed if the sleep deck and mattress are not completely flat.
- Do not load the bed beyond the recommended safe working load of 568.8 lb.
- Do not use the electrical functions of the bed in excess of their designed mode of operation (duty cycle), two minutes on / 18 minutes off.

Important:

As a safety feature, when the actuators of the bed are used too intensively, the bed will go into ‘safe mode’ and shut down. The bed will return to its normal functions after returning to normal temperature. This return to normal temperature may take several hours.

- If the mains cable is in use, ensure that it is not trailing beneath the bed when the bed is moved, as this could result in damage to the cable.
- Always unwind the mains cable fully when it is connected to the mains power outlet.
- Maintenance of the bed must be performed only by qualified personnel.
- When moving the bed, ensure that there is enough space to allow for the movement, in order to avoid damage to the bed, its environment, or causing injury to any person.
- Keep the space around the bed free from any obstructions so that it can be articulated into any of its full range of movements without causing damage to anything or any potential injury to anyone.
- When preparing to move the bed elsewhere, ensure that there is ample space to operate the bed without any obstructions that could create a hazard or cause damage to anything or any potential injury to anyone.
- The electrical and electronic fittings and components used on the bed are not explosion proof and therefore the bed is not suitable for use in hazardous areas, e.g. the bed must not be used in a full oxygen tent.
- Motorised bed mechanisms can cause serious injury.
- Before use of the bed, assess the needs of, and potential risks to, each individual patient.
- Do adopt special safety measures when the bed is to be used by confused patients, e.g. put the bed in a lower position, use safety mats, use protective covers and barriers, enhanced surveillance.
- Movement of the bed could cause damage to soft or improperly installed floor coverings.
- For security of patients, the mains cable should be disconnected and stored safely away from potential harm.
- Only accessories supplied or authorised by Medstrom should be used and then in conjunction with an appropriate risk assessment for patient safety.
- Any serious incident that occurs in connection with the bed should be reported both to the manufacturer and the competent authority of the Country or Member State in which the facility is located.

Avoid Ignition Sources:

- Don't smoke in bed/on the mattress.
- Don't burn candles in the same room as the bed/mattress.
- Don't use matches or lighters in the vicinity of the bed/mattress.
- Don't have electrical equipment in the vicinity of the bed/mattress, e.g. a television over the bed.

- Don't use electric blankets in combination with the bed/mattress.
- Don't have fires and heaters in the vicinity of the bed/mattress.
- Don't place hot items such as hair dryers or heated appliances on the bed/mattress.
- If you use a mobility aid keep it within reach of your bed or device.

5.2. Bed Frame Compatibility



Caution:

Where the mattress used on the bed is not supplied by Medstrom, the compatibility with the bed frame should be assessed independently.

5.3. Contraindications

This medical device is not intended to be used for:

- Patients weighing less than 88.2 lb.
- Patients weighing more than the safe working load (SWL).
- Patients less than 57.5" in height.
- Patients with a Body Mass Index (BMI) of less than 17.
- The bed is not intended for use in ambulances or other mobile medical units.
- The bed is not designed for surgical procedures or as an operating table
- The bed contains materials and components not suitable for MR or certain diagnostic imaging rooms.

If any contraindication(s) exists for an individual patient, carry out and record a risk assessment before determining, as a minimum:

- Is this medical device a suitable and safe choice for the patient?
- Should the patient have access to the wireless handset?
- Are any other safety actions required?

5.4. Potential Risk of Entrapment

The bed has been designed in accordance with IEC 60601 standards to reduce or eliminate the risk of entrapment and/or injury. However, please review the pictorial representations below, which highlight those parts of the bed where small risks remain.



The gap between side panel and sleep deck.



The gap between back rest and side panel.



The gap between the underside of the head/foot board and the floor.



The gap between the underside of the side panel and the floor.

5.5. Installation and Commissioning



Warnings:

Position the bed in such a way that there is always access to the mains power outlet. Do not pull on the mains cable to remove the plug from the mains power outlet.

Before using the bed, ensure that the power requirements of the bed are compatible with the available mains power outlet. The power requirements are provided on the bed reference label and in section “7.2. Technical Specifications”.

New beds are partly dismantled for ease of shipment and require some reassembly prior to use. If the installation is not carried out by Medstrom, proceed as follows:

1. Remove the protective wrapping, securing straps and cable ties.
2. Check that all pin connectors are present and firmly connected to their connections and that the cables are not pinched. Unused connections are sealed to prevent damage from electrostatic discharge.



Warning:

If any electrical assembly or wiring is replaced or repaired, the appropriate electrical safety checks must be made by an authorized individual before returning the bed to use.

3. Thoroughly clean the bed, in accordance with the procedures in section “9. Operational Maintenance”.
4. Position the bed in its use location, ensuring there is enough space to accommodate all bed movements while avoiding damage to the bed or the environment and/or injuring the patient, caregiver or any visitor.
5. Check that the mains cable will not be strained when connected to the mains power outlet.
6. Connect the mains cable, located at the head end of the bed, to the mains power outlet.

Note: For UK use, the cable is provided with a BS1363A, 5A, 3-pin plug.

Note: An equipotential terminal, used to earth the bed, is situated on the frame under the sleep deck.

7. Pair the wireless handset to the control box of the bed, see section “8.1.1. Pairing the Wireless Handset”.
8. Using the wireless handset, check that all functions are operating correctly.
9. Remove mains cable and store in a safe place.

Warning:

When routing cables from other equipment around the bed, check that those cables are not, and cannot be, squeezed, pinched or squashed by any part of the bed.

Caution:

Be aware that:

- You should not change the configuration of the bed without the written permission of the manufacturer.
- Users, patients and their families should be informed of the safety rules to be followed when using the bed.
- The bed comes with a quick user guide (QUG) and training video, these can be accessed via the QR code link on the frame of the bed.

Caution: Electromagnetic Compatibility

Precautions must be taken regarding electromagnetic compatibility (EMC).

- Interference with the bed's electrical/electronic system may be caused by the proximity of communications equipment, mobile phones or other medical equipment such as magnetic resonance imaging (MRI) or active high frequency surgical equipment.
- Caution should be exercised with regard to the use of other electronic equipment in close proximity to the bed. If such usage is unavoidable, the bed and the adjacent electronic equipment should be closely monitored to verify normal operation. If abnormal operation is observed, cease usage of the electronic equipment until an acceptable configuration has been determined and verified.
- The use of cables and accessories other than those specified by Medstrom may negatively affect EMC performance.

Caution: Lithium-ion Battery Pack

The bed is battery powered using a rechargeable long-life lithium-ion battery pack.

When a bleeping sound is heard whilst operating the bed, this indicates a low charge condition. To recharge the battery, connect the bed to the mains power outlet for a minimum of 6 hours. (A full recharge can take up to 10 hours)

6. Transport and Storage



Warning:

This product is not designed for the general transportation of patients, it is designed for patient transportation within the application environments stated in section “3. Intended Use of the Bed”. For general patient transportation use equipment that is specifically designed for that purpose.

- The bed is designed for use in temperatures between 41°F and 104°F.
- The sleep deck must be set flat, the bed height should be set to the low position, all bed movement functions should be locked-out, and the brake must be applied.
- The bed must be shut down as stated in section “6.1. Shutdown Procedure”.
- Use suitable padding materials to protect the bed from the effects of impact and friction.
- During transport in a vehicle, the bed should be restrained to prevent movement.
- Do not use a forklift.
- Beds should not be stacked on top of each other.
- No load should be placed on top of a bed during storage.

6.1. Shutdown Procedure


1. If the mains cable is in use, disconnect it from the mains power outlet and store the cable in a safe place.
2. Lock out all bed movement functions as described in section “8.1. Wireless Handset”.

7. Bed Features

Certified characteristics are safety, electromagnetic compatibility, mechanical safety and fitness for use of the bed and its accessories.

Standard configuration:

- Wireless handset with electric height, back rest, knee rest adjustment and magnetic lock-out.
- One button cardiac chair position and back to flat.
- Electric Trendelenburg/Reverse Trendelenburg.
- Electric CPR.
- Custom height setting.
- Bilateral manual CPR function.
- 2D movement of the back rest by 9° (eliminates patient migration and associated shear and friction on the patient's skin).
- Non removable sleep deck sheets and anti-intrusion plates to prevent access under the bed.
- Four multi-directional double castors x 4” (one anti-static) with front access brake and head end AutoSteer™.
- Long-life battery pack.
- Fixed solid side panel.
- Half-height board secured at the head and foot ends of bed.



Number	Name
1	Back rest
2	Knee rest
3	Leg section
4	Anti-intrusion plates
5	Brake bar

Number	Name
6	Fixed, solid side panels
7	Fixed footboard (half-height)
8	Fixed headboard (half-height)
9	CPR handle

7.1. Dimensions

Description	Dimensions
Height of sleep deck in high position	32.7"
Height of sleep deck in low position	7.5"
Total length	88.6"
Total width	40.6"
Dimension of sleep deck	78" x 35"
Recommended mattress size	78" x 35.4" +/- 0.8"
Half-board height from sleep deck	5.5"
Weight of bed (no mattress or accessories)	304.2 lb
Back rest angle	70°
Leg section angle	20°
Knee rest angle	30°
Trendelenburg angle	+14°
Reverse Trendelenburg angle	- 14°
Patient height	Minimum 57.5" / maximum 72.8"
Minimum patient weight	88.2 lb



Warning:

Incompatible mattresses can create hazards. Specialist mattresses require a local risk assessment.

Many specialty mattresses feature immersion therapy which allows the patient to sink into the surface of the mattress, thus reducing the risk of patient falls.

7.2. Technical Specifications

Description	Dimensions
Mechanical	
Safe working load	568.8 lb (425.5 lb patient, 143 lb for a mattress plus bedding and accessories)
Castors	Four multi-directional double castors x 4" (one anti-static)
Brake & steer	Brake at foot end, steer at head end
Electrical	
Power Input	Voltage 100-240 V ac 50/60 Hz
	Current 3.9 A Max
Mode of operation (duty cycle)	Non-continuous: 2 min. on/18 min. off
Output voltage of the transformer	24 to 26 VCC
Amperage	5 A max per channel, 5 A max total on secondary
Electric power	0.5 W idle, 200 W max full load
Electronic protection	Individual in case of overload. Durability of non-resettable thermal circuit breaker relays.
Fuses	Electronic overcharge protection
Wireless handset	Linak type HB200
Power supply unit	Linak type CO61-BT-Mk2
Power supply wire	3x1.5mm ² PVC sleeve– interchangeable
Variable height motor	Linak type LA40 – 6000N
Back rest motor	Linak type LA27 – 3500N
Knee rest motor	Linak type LA27 – 3000N
Protective grounding class	Class I
Class of protection against dust and water	IPX4 for wireless handset HB200
	IPX6 for control box
Classification of the device	Type B according to IEC 60601
Battery	BA22 Lithium-ion
Environmental Conditions of Use and Storage	
Use temperature	From 41°F to 104°F
Storage temperature	From 23°F to 122°F
Thermal cut-out reset time	Up to 5 hours
Humidity	From 20% to 90% - to 30° without condensation
Atmospheric pressure	700 to 1060 hPa
Maximum temperature of parts	139.6°F (for electrical boxes)
Sound level	41 dB

8. Instructions for Use



Warnings:

For your safety, before using the bed, it is essential to verify that:

- The mattress you intend using is compatible with the bed.
- The recommended mattress dimensions are 78" x 35.4", +/- 0.8".
- The recommended mattress depth is 5.5". Maximum depth of 9.8" when using a specialist mattress.

Ensure the brake and steer are applied properly in the locked position when the bed is stationary.

Mattress and side panel combinations that are incompatible with the dimensions that are provided in this document can cause an entrapment hazard.



Warning:

Keep clear of the bed when it is being operated. Severe injury can result from crushing by moving parts.

Be aware that bed sections can be operated under battery power even when the mains cable is disconnected from a mains power outlet.

8.1. Wireless Handset

8.1.1. Pairing the Wireless Handset

In order to use the wireless handset, it first needs to be paired with the bed using Bluetooth[®].


1. Stand within 1 metre of the head end of the bed.
2. Enter direct pairing mode by holding the magnetic key over the O symbol and pressing the "bed down" button at the same time and hold for three seconds, until a beep is heard.
3. Release the "bed down" button maintaining the magnet in place and move towards the control box – the beeping frequency will speed up.
4. Confirm pairing by pressing the "bed down" button again.
5. Test the pairing is successful by pressing any of the bed functions to check movement wirelessly.



 Warning:

Keep clear of the bed when it is being operated. Severe injury can result from crushing by moving parts.
Be aware that bed sections can be operated under battery power even when the mains cable is disconnected from a mains power outlet.

8.1.2. Wireless Handset



Control	Action
Electric CPR buttons	Press and hold the back rest up and down buttons simultaneously. This will flatten the sleep deck and set the bed to a height suitable to perform CPR.
Back rest buttons	Press and hold the up/down button to raise or lower the back rest to the desired angle.
Bed height buttons	Press and hold the up/down button to raise or lower the sleep deck.
Knee rest buttons	Press and hold the up/down button to raise or lower the knee rest to the desired angle.
Cardiac chair/back to flat buttons	Cardiac chair: - Press and hold the up cardiac chair button. The back rest will raise to 45° and the knee rest will raise to 30°. The bed will then go into the reverse Trendelenburg position. Back to flat: - To return the bed to the flat position, press and hold the down cardiac chair (back to flat) button.

Continued...

Control	Action
Custom bed exit height setting	Press and hold the cardiac chair button and back to flat button simultaneously for 10 seconds to set the correct bed exit height for an individual patient (see section “8.5.6. Storing a Custom Bed Exit Height in the Memory”).
Trendelenburg / Reverse Trendelenburg buttons	Press and hold the Trendelenburg button. This will flatten the sleep deck, pause it for two seconds, then tilt the sleep deck head downwards utilising the height actuators. Press and hold the reverse Trendelenburg button. This will flatten the sleep deck, pause it for two seconds, then tilt the sleep deck head upwards utilising the height actuators.
Function lock-out	To lock out the handset, hold the magnet over the circle symbol just above the Medstrom logo and press the reverse Trendelenburg button until the battery light on the left illuminates (after two or three seconds). The handset is then locked (light will illuminate when pressing any function and that function will not move). To unlock, use the magnet again and press the Trendelenburg button until the battery light goes out (after two or three seconds). The handset is then unlocked.

Note: It is necessary to keep pressing the function button for the entire duration of the desired movement. An exception to this is when lowering the bed height to its minimum, the bed goes down and reaches a height of 9.4” from the floor. At this point the bed pauses to indicate a risk of entrapment if an object or bodily part is located between the bed and the floor. To continue lowering to the lowest position, release the button, then press it again. This will lower the bed to the ultra-low position. An audible signal will sound to alert the user of the reduction of the space between the floor and the bed.

8.2. Recharging the Battery

A fully charged battery will provide sufficient power to operate one full bed movement each day for 28 days.

- 1. Locate the charging point under the backrest, patient right-side. (Label indicates where)
- 2. Connect mains cable to the charging point and plug the 3-pin plug into the mains power socket.



Note: For a very low or empty battery allow up to 10 hours to fully charge.

8.3. Variable Bed Height Function

When the down bed height button is pressed and held, the bed will lower and stop at a pre-set height of 9.4”, before it reaches its lowest position. The bed will sound an alarm after it reaches the 9.4” position.
If a custom height has been programmed, the bed will stop at that custom height and at 9.4”.
In its lowest position, the sleep deck is 7.5” from the floor. The bed can be positioned at any height between the highest and lowest positions.
Maximum care height for the bed is 31.9”.





The brake will engage automatically when the bed reaches its lowest position. The bed must be raised by 2" from the lowest height (7.5" + 2") in order to access the brake and allow the bed to be moved. The lowest height position of the bed should be risk-assessed for patients at high risk of falling.

It is recommended that periodically the bed is placed into either the highest or lowest position to maintain the precision of the height adjustment software.

Operation:

The height of the bed can be adjusted using the wireless handset. Press the up button to go up or the down button to go down to the desired height.



Warnings:

For patients who may be at risk of injury from falls, carry out a risk assessment to determine if leaving the bed at the low height may reduce the risk of injury.


Under certain circumstances it may be necessary to prevent patient operation of the bed functions, this can be done using the selective lock out function on the patient handset.

Hoists and under-bed tables cannot fit under the bed when it is in its ultra-low position.

8.4. Sleep Deck Articulation Functions

The sleep deck can be articulated using the following functions:

- Back rest function.
- Manual CPR function.
- Electric CPR function.
- Knee rest and leg sections function.
- Cardiac chair function.
- Storing a custom bed exit height in the memory.
- Trendelenburg / reverse Trendelenburg function.



Articulated sections of the sleep deck

8.4.1. Back Rest Function

This allows angle adjustment of the back rest.

The back rest will pause at the 30° and 45° angles when it is raised or lowered. This eliminates the need for angle indicators.

The 2D back rest design eliminates patient migration, particularly when the cardiac chair button is used. This gives better comfort to the patient and reduces the risk of shear and friction.

The back rest angle is adjustable from 0 to 70°.

Electric adjustment of the back rest only occurs with the use of the wireless handset.

Operation:

- Press the back rest up button to raise the back rest, or the back rest down button to lower the back rest.



Caution:

This bed is equipped with both manual and electric CPR functions. The manual CPR function uses a damper to assist with controlling the back rest movement. The electric CPR function is operated by simultaneously pressing the two (up/down) back rest buttons.

8.4.2. Manual CPR Function

The manual CPR function is a quick way of setting the back rest horizontal in readiness to perform CPR.

Operation:

The following operation should only be performed if the bed is occupied. The electric back rest function will operate normally.

- Pull and hold either of the orange CPR handles and lower the back rest until flat.
- Release the CPR handle.



Bilateral CPR handles

Caution:

This function should not be used as a simple back rest tilt adjustment. Improper use of the manual CPR release handles may degrade and damage the mechanism.

8.4.3. Electric CPR Function

The electric CPR function flattens the sleep deck within seconds and sets the bed to a height suitable to perform CPR.

Operation:

- To operate the electric CPR function using the wireless handset, press and hold the up and down back rest buttons simultaneously.
- This will set the sleep deck to the flat position and to a height appropriate for CPR.



8.4.4. Knee Rest and Leg Sections Function

This function allows adjustment of the knee rest and leg sections.

Operation:

- Press the up knee rest button to raise the knee rest, or the down knee rest button to lower the knee rest.



8.4.5. Cardiac Chair Function

The cardiac chair function allows the patient to be placed in a full chair position. It combines movement of the back rest and the knee rest with the reverse Trendelenburg position in a continuous motion, using only one button.



Operation:

- Press and hold the up cardiac chair button using the wireless handset. The back rest will raise to 45° and the knee rest will raise to 30°. The bed will then go into the reverse Trendelenburg position.
- To revert the bed to flat, press and hold the down cardiac chair (back to flat) button.

8.4.6. Storing a Custom Bed Exit Height in the Memory

This function allows the 'bed exit position' to be set for each patient. This promotes safe patient mobilisation. The exit position ensures that the bed is at the correct height for the patient to get their feet flat on the floor prior to raising themselves to exit the bed. This greatly reduces the risk of falls during mobilisation.

It is recommended that a custom exit position is set for every patient when they are first placed on the bed.

Operation:

- Using the wireless handset, put the bed at the height where the patient's feet are flat on the floor and their hips and knees are at an approximately 90° angle.
- Then, simultaneously press the up cardiac chair button and the down cardiac chair (back to flat) button for 10 seconds.
- A beep will sound when the bed height position is registered and stored in the memory. The bed will now stop at this position when raising or lowering the bed until the process is repeated for the next patient.



8.4.7. Trendelenburg/Reverse Trendelenburg Function

This function should be used with the sleep deck in the flat position. Tilting the sleep deck without first setting it back to flat will provoke a downward sliding of the patient.

Operation:

The Trendelenburg (feet up) function allows you to tilt the whole bed backward and down.

- Using the wireless handset press and hold the down cardiac chair (back to flat) button to set the sleep deck to the flat position.
- Then press and hold the right-hand up Trendelenburg/reverse Trendelenburg button until the bed reaches the desired feet up angle.
- To return the sleep deck to the flat position, press and hold the left-hand up Trendelenburg/reverse Trendelenburg button. The function will pause to indicate when the sleep deck is back to horizontal. Release the button.
- Continuing to hold the left-hand up Trendelenburg/reverse Trendelenburg button would allow the reverse Trendelenburg (feet down) movement to begin.

The reverse Trendelenburg (feet down) function allows you to tilt the whole bed forward and down.

- Using the wireless handset, press and hold the down cardiac chair (back to flat) button to set the sleep deck to the flat position.
- Then press and hold the left-hand up Trendelenburg/reverse Trendelenburg button until the bed reaches the desired feet down angle.
- To return the sleep deck to the flat position, press and hold the right-hand up Trendelenburg / reverse Trendelenburg button. The function will pause to indicate when the sleep deck is back to horizontal. Release the button.
- Continuing to hold the right-hand up Trendelenburg / reverse Trendelenburg button would allow the Trendelenburg (feet up) movement to begin.




Caution:

The Trendelenburg / reverse Trendelenburg positions are related to clinical indications and should be used only upon advice of a licensed healthcare professional. Only use the Trendelenburg / reverse Trendelenburg function when the brake is locked. Ensure there is sufficient head space around the bed to accommodate the manoeuvre.

8.4.8. Lockout of Functions

To lock out the handset, hold the magnet over the magnet symbol, just above the Medstrom logo, and press the reverse Trendelenburg button until the battery light on the left illuminates (after two or three seconds). The handset is then locked. The light will illuminate when pressing any function and that function will not move.



To unlock, use the magnet again and press the Trendelenburg button until the battery light goes out (after two or three seconds). The handset is then unlocked.



8.5. Braking System

The brake is located at the foot end of the bed.

The brake should be applied (locked/pedal in the down position) at all times that the bed is stationary.

The bed has been designed such that the brake is applied automatically when the bed height is set to its lowest position of 7.5". Whilst in this position it is not possible to release the brake. In order to access the brake, the bed must be raised by at least 2" above the lowest height setting (7.5" + 2").

Operation:

- To apply the brake, using one foot, press down the brake bar at the foot end of the bed.
- To release the brake, using one foot, insert the toes under the brake bar and lift up (AutoSteer™ mode).



Brake released – brake bar in up position

Caution:

The brake function of the bed satisfies the requirements of the basic safety and essential performance standard IEC 60601-2-52. This requires the brake system to be certified to hold the bed stationary when it is at an angle of 6° relative to the horizontal plane on a concrete floor covered with a vinyl coating of 1/16" to 1/8" thick, in normal use.

Caution:

Be mindful of possible trapping of a foot when setting the bed height to its lowest position as the brake bar is automatically applied in the lowest position.

**Warning:**

Maintenance of the brake should only be carried out by an authorised technician.

**Caution:**

If the bed does not have to be moved, it is strongly advised to apply the brake to avoid potential falls when a patient is accessing or exiting the bed.

Before stopping the bed in its final position, ensure that nothing will prevent normal movement of the bed; up, down or sideways.

8.6. Steering System - AutoSteer™

As standard the bed is equipped with eight multi-directional castors, in four castor assemblies.

The steering lock acts on the front set of castors on the castor assemblies at the head end of the bed. The steering lock is applied or released using the AutoSteer™ pedal at the head end of the bed.

- When the steering lock is applied, AutoSteer™ pedal down, the front castors are locked in the straight-ahead position. This allows for easier control of bed movement and prevents lateral (sideways) movement of the bed (AutoSteer™ mode).
- When the steering lock is released, AutoSteer™ pedal up, the front castors are released. This allows the bed to be moved in any direction, including lateral (sideways) movement (FreeMove™ mode).

The brake and steering should be applied (locked/pedals in the down position) at all times that the bed is stationary.

Operation:

- To apply the steering lock, using one foot, press down the AutoSteer™ pedal at the head end of the bed.
- To release the steering lock, using one foot, insert the toes under the AutoSteer™ pedal and lift up.



Steering lock applied – AutoSteer™ pedal in down position

8.7. Moving the Bed

It is recommended that when moving the bed it should be pushed and steered from the foot end of the bed.

When travelling in a straight line along corridors or when leaving the bed unattended, the steering lock should be applied, by putting the AutoSteer™ pedal in the down (locked) position. This will allow controlled movement of the bed with reduced effort, significantly reducing the force required to move the bed, and helps steer the bed around tight corners, with or without the patient being in the bed (AutoSteer™ mode).

When lateral (sideways) movement of the bed is required, the steering lock should be released, by putting the AutoSteer™ pedal in the up (released) position. This will allow movement of the bed in any direction (FreeMove™ mode).

To locate the bed into a bed space, make sure steering lock is applied, AutoSteer™ pedal in the down (locked) position, and push the bed into place in a straight line. When the bed is in position, apply the brake. Ensure the brake bar is in the down (locked) position.

The brake and steering should be applied (locked/pedal in the down position) at all times that the bed is stationary.

Operation:

- For small movements of the bed, if the mains cable is still connected to the mains power outlet, check that the mains cable is not being pulled tight and being stressed.
- For larger movements, where the mains cable is still connected to the mains power outlet, disconnect the mains cable from the mains power outlet, and stow it safely to avoid it trailing on the floor.
- Check that no device or accessory is connected to the mains power outlet before moving the bed.
- Release the brake, with brake bar in the up (released) position, and leave the steering lock applied, AutoSteer™ pedal in the down (locked) position.
- Manoeuvre the bed by holding either the headboard or footboard, ideally the footboard.
- When locating the bed in a bed space or patient room, check that no obstacle will prevent normal movement of the bed; up, down or sideways.



Caution:


It is advised to position the bed at least 13.8" from the wall during installation. This will ensure that if the Trendelenburg/reverse Trendelenburg function is used, the bed will not touch the wall.

9. Accessories:

- Fixed Castor Kit

9.1. Mental Health Fixed Castor Kit

A Mental Health Fixed Castor Kit is available to replace the standard castors with two fixed 4" Castors at the head end and two Feet at the foot end which can be fixed to the floor.



Mental Health Fixed Castor Kit

10. Operational Maintenance



Warning:

The bed is not designed for automated cleaning - only manual cleaning should be undertaken.



Caution:

Using detergent/disinfectant products does involve some risk. Always wear personal protective equipment (PPE), follow the instructions and avoid mixing different products.

10.1. Cleaning Guidelines

The bed should be cleaned and disinfected weekly and between patients, following these guidelines:

- Unplug the bed (if the mains cable is connected), and lock out all of its electric functions using the wireless handset (see section "8.5.8. Lockout of Functions").
- Ensure that the mains cable plug does not come into contact with any liquid.
- Check that electrical connections are secure.
- Check that electrical components do not show any signs of wear likely to allow the penetration of liquid.
- Do not wash the bed in a wash tunnel, using a water jet or with any sort of pressure equipment.
- Avoid saturating areas containing electrical connections or electrical components. Once cleaned, dry areas carefully to avoid any risk of moisture remaining, in particular in the area of electrical connections.



Warning:

If there is any doubt or concern about the penetration of liquid into electrical connections or electrical components, it is strongly recommended that the bed should not be reconnected to the mains power outlet and the maintenance department advised.

- Clean surfaces with a soft cloth, warm water and a mild detergent/disinfectant solution. Ideally use the cleaning/disinfection solution recommended by the care facility, following the dilution instructions carefully.
- The use of abrasive products or materials is forbidden.
- Products such as petrol, ketonic solvents, concentrated alkali or acid products, chlorine-based solvents or dissolvent risk causing permanent damage to the surfaces and should not be used.
- Non-ionic detergents diluted to 5%, bleach diluted to 10%, and ammonia-based products diluted to 5% may be used.
- Stains left by coloured substances such as eosin, betadine, etc., as well as food products, should be cleaned as quickly as possible to avoid the risk of contamination to the patient.
- For stubborn stains, the use of pure disinfection products is possible locally, as long as the necessary precautions are taken.
- Areas showing traces of cuts or deep scratches should be repaired to avoid the risk of infiltration and deterioration of the protective surface.

10.2. Sources of Dirt and Contamination

Sources of dirt and contamination of a medical bed in use are:

- Fabric dust coming from the bedding.
- Food residues.
- Ink.
- Antiseptic liquids.
- Vomit, urine, excrement, blood, etc.

It is necessary to take three types of cleaning/disinfection or bio-cleaning of the bed into consideration:

1. Daily bio-cleaning of the areas of high risk.
2. Bio-cleaning at the start and on the transfer of the patient, and, as a minimum, every month for high and moderate risk areas.
3. Full bio-cleaning of the bed after the departure of a patient presenting risk of infection, and, as a minimum, every two months.

10.3. Recommended Cleaning Procedures

Healthcare establishments are divided into sectors in terms of the risk of infection.

Frequency and methods of cleaning and disinfection are adapted according to the evaluated risk.

The level of infection risk in one room with a patient being cared for is moderate. It is much higher if the patient presents an identifiable infection risk.

Cleaning or disinfection of a medical bed, although not in direct contact with the patient and in particular damaged areas, should be performed regularly. Certain areas of the bed are more subject to contamination given the frequency of contact with caregivers' and patients' hands. The bed can be broken down into three areas:

High Contamination Level Areas

- Handle areas of the headboard and footboard.
- Fixed, solid side panels.
- Wireless handset.

Moderate Contamination Level Areas

- Bed headboard and footboard.
- Upper surfaces of the bed.
- Brake bar and AutoSteer[™] pedal.
- Anti-intrusion panels.

Low Contamination Level Areas

- Metal structure of the bed.
- Underside surfaces of the bed.
- Electric actuators.
- Wheels.

10.3.1. Daily Bio-Cleaning Procedure

This procedure can be carried out with the patient still in the bed. The objective being to ensure good hygiene of the parts regularly in contact with caregivers' and/or patients' hands.

- Carefully clean the handle areas of the headboard and footboard, or fixed, solid side panels and the wireless handset.
- Eliminate all traces of dirt apparent on the other parts of the bed.

Important:

Clean the side rail in the upright position.

10.3.2. Monthly Bio-Cleaning Procedure or Bio-Cleaning Procedure on the Departure of a Patient

This procedure is carried out without a patient in the bed. The objective being to disinfect all parts of the bed that regularly come into contact with hands, plus areas dirtied by liquid deposits, secretions, dust and food residues, etc.

- Use locally approved cleaning materials and detergent/disinfectant solutions.
- Release the brake and move the bed away from the wall to give access all-round the bed; then reapply the brake.
- Put the bed in the flat position and set to a suitable working height.
- If the bed is plugged in, remove the mains cable from the mains power outlet and lock out all electrical functions using the wireless handset.

Cleaning the Foot End of the Bed

- Clean the footboard.
- Lift the foot end of the mattress and fold over towards the head end of the bed. Clean the upper surfaces of the bed, including the anti-intrusion plates, and the inner side of the fixed, solid side panels and then the underside of the mattress.
- Fold the mattress back into place. Clean the outer side of the fixed, solid side panels.

Cleaning the Head End of the Bed

- Clean the headboard.
- Clean the wireless handset and both orange CPR release handles.
- Lift the head end of the mattress and fold over towards the foot end of the bed. Clean the upper surfaces of the bed, including the anti-intrusion plates, and the inner side of the fixed, solid side panels and then the underside of the mattress.
- Fold the mattress back into place. Clean the outer side of the fixed, solid side panels.
- For beds with Autosteer™, clean the AutoSteer™ pedal.
- Release the brake and re-position the bed against the wall; then reapply the brake.
- Clean the brake bar.
- Reconnect the bed to the mains power outlet (if required) and unlock the function lock outs using the wireless handset (see section "8.5.8. Lockout of Functions").

10.3.3. Full Bio-Cleaning Procedure for the Bed

This procedure is carried out without the patient in the bed. The objective is to disinfect the whole bed after it has been occupied by a contaminated patient or periodically every two months.

This operation should also be carried out before first using the bed.

- Use locally approved cleaning materials and detergent/disinfectant solution.
- Release the brake and move the bed away from the wall to give access all-round the bed; then reapply the brake.
- Remove the mattress.



Warning:

The central part protects the power supply unit, its connections and two sensors.

- Adjust the bed base to a suitable working height and raise all the hinged sections using the wireless handset.
- If the bed is plugged in, remove the mains cable from the mains power outlet and lock out all electrical functions using the wireless handset (see section “8.5.8. Lockout of Functions”).

Cleaning the Upper Part of the Bed

- Clean the footboard and headboard, the wireless handset and the sleep deck and anti-intrusion plates all around the deck.

Cleaning the Lower Part of the Bed

- Clean the bed chassis, the main frame and both orange CPR handles.
- Working down the bed, clean the bed legs, the elevation arms and the castors.
- For beds with Autosteer™, clean the AutoSteer™ pedal.
- Release the brake and re-position the bed against the wall. Reapply the brake.
- Clean the brake bar.
- Reconnect the bed to the mains power outlet, remove the electrical function lock outs, and set the bed to the flat position using the wireless handset.
- Clean the mattress and replace it on the bed.

11. Electromagnetic Compatibility (EMC):

**Warning:**
Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12” to any part of the Solo MH bed, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The purpose of this testing is to ensure the Solo MH bed is not likely to adversely affect the normal operation of other such equipment and that other such equipment is not likely to adversely affect the normal operation of Solo MH bed.

**Warning:**
Normal operation of the Solo MH bed can be affected by using a mobile or a microwave oven. HF surgical equipment, magnetic resonance imaging or other radio radiant equipment near this product may cause malfunction or lead to loss of essential performance.

Table 2-1: Guidance and Manufacturer’s Declaration— Electromagnetic Emissions (IEC 60601-1-2)		
The Solo MH bed is intended for use in the electromagnetic environment specified below. The customer or the user of the Solo MH bed should ensure that it is used in such an environment.		
Emissions Test	Compliance Level IEC 60601-1-2	Electromagnetic Environment Guidance
Radiated Emissions FCC Part 15:2016 CISPR 11:2015+A2:2019	Class A	The Solo MH bed uses RF energy only for receiving control functionality. Testing completed shows that RF emissions are low and do not cause any interference with nearby electronic equipment.
Conducted Emissions (Minimum) FCC Part 15:2016 CISPR 11:2015+A2:2019	Class A	The Solo MH bed is suitable for use in all professional healthcare establishments directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Conducted Emissions (Maximum) FCC Part 15:2016 CISPR 11:2015+A2:2019		

Table 2-2: Guidance and Manufacturer's Declaration — Electromagnetic Immunity (IEC 60601-1-2)

The Solo MH bed is intended for use in the electromagnetic environment specified below.

The customer or the use of the Solo MH bed should ensure that it is used in such environment. The impact on electromagnetic interference from unique medical emitters such as electrocautery devices, electrosurgical units, and diathermy equipment has not been determined. These EM emitters should be avoided.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance	
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	Compliance level met or exceeded	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.	
Electrical Fast Transient/Burst IEC 61000-4-4	± 2 kV for power supply lines	Compliance level met or exceeded	Mains power quality should be that of a typical commercial and/or hospital environment.	
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to earth	Compliance level met or exceeded	Mains power quality should be that of a typical commercial and/or hospital environment.	
Voltage Dips, Short Interruptions, and Voltage Variations on Power Supply. IEC 61000-4-11 (see Note 1)	< 5 % UT (> 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (95% dip in UT) for 5 seconds	Compliance level met or exceeded	Mains power quality should be that of a typical commercial and/or hospital environment. If the user of the Solo MH bed requires continued operation during power mains interruptions, it is recommended that the Solo MH bed be powered from an uninterruptible power supply or a battery.	
Power Frequency Magnetic Field IEC 61000-4-8	30 A/m (50/60 Hz)	Compliance level met or exceeded	Professional Healthcare Levels	
Proximity Magnetic Fields Immunity IEC 61000-4-39	30 kHz CW 8A/m 135.2 kHz PM 65 A/m 13.56 MHz PM 7.5 A/m	Compliance level met or exceeded	Professional Healthcare Levels	
Conducted RF IEC 61000-4-6	150 kHz – 80 MHz 6Vrms outside ISM bands* 150 kHz – 80 MHz 6Vrms in ISM bands*. (V1)	Compliance level met or exceeded	Portable and mobile RF communications equipment, including cables, should be used no closer to any part of the Solo MH bed than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d=2\times\sqrt{P}$	
Radiated RF IEC 61000-4-3	80 MHz – 2.7 GHz 3 V/m (E1) 380–390 MHz 27 V/m 430–470 MHz 28 V/m 704–787 MHz 9 V/m 800–960 MHz 28 V/m	Compliance level met or exceeded	$d=4\times\sqrt{P}$ 80 MHz to 800 MHz	$d=7.7\times\sqrt{P}$ 800 MHz to 2.7GHz
			Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.	

Table 2-2: Guidance and Manufacturer's Declaration—Electromagnetic Immunity (IEC 60601-1-2)

The Solo MH bed is intended for use in the electromagnetic environment specified below.

The customer or the use of the Solo MH bed should ensure that it is used in such environment. The impact on electromagnetic interference from unique medical emitters such as electrocautery devices, electrosurgical units, and diathermy equipment has not been determined. These EM emitters should be avoided.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
	1700–1990 MHz 28 V/m 2400–2570 MHz 28 V/m 5100–5800 MHz 9 V/m		 Interference may occur in the vicinity of equipment marked with this symbol: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12" to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
NOTE 1:	UT is the AC mains voltage before application of the test level.		
NOTE 2:	At 80 MHz and 800 MHz, the higher frequency range applies.		
NOTE 3:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.		

* The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Solo MH bed is used exceeds the applicable RF compliance level above, the Solo MH bed should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Solo MH bed. Over the frequency range 150 kHz to 80 MHz, field strength should be less than 3 V/m.

12. Warranty, Servicing and Maintenance

Medstrom reserves the right to make changes to the design, characteristics and models without prior notice.

The only warranty Medstrom makes is the express written warranty extended on the sale or rental of its products.

Warranty may lose its validity in the following cases:

- Disassembly of mechanical or electrical parts of the bed without consultation and approval of the manufacturer.
- Replacement parts are not supplied by the manufacturer.
- Work on the electrical parts or gas cylinders by un-approved personnel.
- Degradation of coatings or materials from shock, friction and scratches.
- Abnormal use of bed not following the precautions and recommendations, including SWL and maximum/minimum patient weights, resulting in degradation of the bed or its environment.
- Medstrom do not approve disinfectants. Disinfectants are locally approved by each establishment.
- Use of cleaning products without respecting appropriate dilution proportions.
- Washing with a water pressure jet.
- Crushing or cutting the mains cable or the control cables.
- Intensive use of electric bed functions beyond the recommended service factor.
- Use of a bed that indicates a mechanical or electrical malfunction.
- Stacking of beds during storage.
- The manufacturer, assembler or installer cannot be held responsible for the safety, reliability or the characteristics of the device if it is not used in a safe electrical environment.
- Accommodating the bed complies with the relevant recommendations, and the device is used in accordance with the instructions for use.



Warnings:

- It is prohibited to modify the bed.
- Use only spare parts approved by the manufacturer.
- Any item damaged must be replaced before using the bed.
- Disposal of the product is necessary if the essential requirements are no longer met, especially when the product no longer has its original features.
- On disposal, product must be rendered unusable.
- Please observe the environmental regulations of the country in which the product is being disposed of.
- The manufacturer specifies that the mattress is not an element of the bed and is not part of this manual.

13. Service Checklist Form

General Condition			
	Pass	Fail	Comment
Castors: Rotate through 360°, check fixings are secure.			
Brake mechanism: Lock in position, condition of brake castors.			
Steer mechanism: Lock in position.			
Side panels: Free from damage and scratches. Check fastenings are secure.			
Frame: Free from distortion, cracks in welds, fastening points secure, spacers.			
Sleep deck: Securing points in place, free from damage, check fastenings are secure.			
Handset: Paired correctly, all functions working, free from damage.			
Headboard/footboard: Securely locked in position, free from damage.			
Labels: All labels applied, i.e. serial number, SWL, model, asset label, brake/steer.			
Electrics: Check cables for wear and tear, damage and secure connections (bed must not be connected to mains power during these checks)			
Functionality			
Raise and lower: Operate to set positions.			
Back rest section: Operate to set positions, pause at 30° and 45° angles.			
Knee rest: Operate to set positions.			
Reverse Trendelenburg/Trendelenburg: Operate to set positions.			
Lateral bed movement: Bed can be moved with ease in lateral direction (FreeMove™).			
Castors: Braking - remains in position when applied.			
Castors: Steer - travels in straight line when steering pedal is engaged (AutoSteer™).			
Features			
Ultra-low position: Operational beyond standard low position.			
CPR: Operational via the handsets, puts bed at working height and in a flat position.			
Manual CPR: Releases back rest section to flat position once handle is pulled.			
Custom height: Sets custom height to designated position from the wireless handsets.			
Bed reset (calibration): Recalibrates control box with actuators operated from handsets (raise bed up to highest, then lowest position. Same process with back rest and knee rests).			
Date:			
Signed:			

14. Manufacturer's Details

UK:



Medstrom Ltd
2 Cygnus Court, Beverley Road,
Pegasus Business Park, Castle Donington,
DE74 2SA
United Kingdom
Tel: +44 (0)345 371 1717
info@medstrom.co.uk

US Distributor:



Pineapple Contracts Inc.
7015 19 Mile Rd, Sterling Heights
MI, 48314
USA
Tel: +1 (800) 496 9324
gspr@pineapplecontracts.com

For further information relating to technical characteristics, maintenance or after sales service, please consult the bed's technical manual available on request.

Product Numbers:

Solo MH:

- Solo MH Side Panel 10P Configuration SOLOMHFGPNL10P-B-USP

GTIN Numbers:

Product	Region	GTIN
Solo Mental Health	US	05060467210829

This page is intentionally blank.

Pineapple

manufactured by **medstrom⁺**

Medstrom Ltd, 2 Cygnus Court, Beverley Road, Pegasus Business Park, Castle Donington, Derby, DE74 2SA

PRO00096_US, Revision 1.0
CRN-0000338

