

# Horizon 10 Adjustable Height Cradle

## Technical Specification

- Cradle made of shatterproof, transparent plastic with aeration holes and padded mattress
- Epoxy coated steel trolley with 4x75mm revolving castors, two with individual diagonal brakes
- Padded mattress made of washable plastic fabric
- Hydraulic adjustable height
- A safety ring on top allows the cradle tilting up to 12°
- **Overall dimensions:** 770mm x 500mm x 1090-1305mm
- **Height from floor to mattress:** 820-1035mm
- **Cradle dimensions:** 770mm x 480mm x 195-240mm
- **Mattress dimensions:** 610mm x 370mm x 50mm
- **Safe Working Load:** 10kg
- **Total weight:** 16.2kg

