

Horizon 10 Fixed Technical Specification



SPECIALIST PRODUCTS • UNIQUE SERVICES

Improving outcomes in the patient environment

- Cradle of neutral transparent, shatterproof plastic with aeration holes and padded mattress
- Frame made of steel tube painted with anti-scratch epoxy powders in RAL 9006, metal grey colour
- Trolley with 4 x 75mm revolving castors, two with individual brake diagonal
- Padded mattress made of a washing fabric. The mattress is hypoallergenic, antibacterial and is resistant to staining
- **Overall dimensions:** 770mm x 480mm x 890mm
- **Height from floor to mattress:** 660mm
- **Trolley dimensions:** 735mm x 440mm x 745mm
- **Cradle dimensions:** 770mm x 480mm x 210/270mm
- **Mattress dimensions:** 610mm x 370mm x 50mm
- **Safe working load:** 10kg
- **Weight:** 12kg

