

# Technical Data:

## Bari10A

**medstrom<sup>+</sup>**

Improved Patient Outcomes



Bed Weight	330kg
Safe Working Load	477kg (75st)
Minimum Height	48cm
Maximum Height	86cm
Overall Bed Dimensions	Minimum 99cm (W) x 229cm (L) / Maximum 132cm (W) x 249cm (L)
Sleep Deck Dimensions	Minimum 91cm (W) x 203cm (L) / Maximum 122cm (W) x 223cm (L)
Back Rest Angle	65°
Thigh Section Angle	35°
Calf Section Angle	10°
Trend/Reverse Trend Angle	+/- 15°
Castors	Heavy duty steer and brake castors
CPR	Electronic and manual CPR controls
Side Rails	Simple folding split rail design
Warranty	1 year

SM1146, Rev2, Feb2025