

Technical Data: P.R.O. Matt Plus Range



Minimum patient weight	No minimum patient weight.
Maximum patient weight	227kg (35.7st)
Mattress width (Std)	85cm
Mattress width (Wide)	107cm
Mattress length	203cm
Mattress depth	17cm
Mattress weight (Std)	11kg
Mattress weight (Wide)	14kg
Control unit dimensions	8cm (w) x 26cm (h) x 13cm (d)
Control unit weight	1.5kg
Comfort settings	3
Modes of operation	Continuous low pressure, alternating (1-in-2 cell), dual modea and auto firm.
Cycle times	5, 10, 15 or 20 minutes.
Lock-out	Yes, automatic after 30 seconds of inactivity.
Alarms	Low pressure and power fail (visual and audible).