


Technical Data: Tilt-in-Space Chair

medstrom⁺

Improved Patient Outcomes



Seat Width	28" / 71cm
Overall Width	36" / 92cm
Seat Depth	24" / 60cm
Overall Length	75" / 190cm (with back rest down & leg rest up)
Seat Height	22" / 55cm (normal position)
Overall Height	52" / 132cm
Seat Back Height	29" / 73cm
Arm Height	6" / 16cm
Overall SWL	55st / 350kg
Leg Rest Length	17" / 44cm (not extended), 22" / 55cm (extended)
Castor Diameter	6" / 15cm (twin wheel castors with brake & steer options)
Seat Composition	Integrated gel seat cushion with lumbar & backrest air cell support.
Powered Options	Back rest, leg rest, chair tilt, chair height



Scan the QR code to access the clinically focused training video

SM609, Rev3, Feb 2021