

Technical Data:

Riser Recliner PRO

medstrom⁺

Improved Patient Outcomes



Seat Width	20" / 51cm	24" / 61cm	28" / 71cm	32" / 81cm	34" / 86cm	38" / 96cm
Overall Width	30" / 77cm	34" / 87cm	38" / 97cm	42" / 107cm	44" / 112cm	48" / 122cm
Seat Depth	20" / 51cm					
Overall Depth	36" / 75cm					
Seat Height	15" / 40cm - 27" / 70cm					
Overall Height	45" / 114cm					
Seat Back Height	28" / 71cm					
Arm Height	7" / 18cm					
Overall SWL	198kg / 31st	198kg / 31st	390kg / 61st	390kg / 61st	444kg / 70st	444kg / 70st
Braking System	Electronic braking - safe braking before any chair articulation delivered					
Hoist Clearance	Chair height can be raised to allow for a minimum clearance of 2" / 5cm					
Upholstery Composition	Multi-stretch, durable outer vinyl with Shield plus for infection prevention					
Powered Options	Chair lift, chair tilt, independent back rest, independent leg rest (with extension)					
Back Rest	Waterfall back rest to reduce shear & friction (with removable lumbar cushion)					

SM1132, Rev2, CRN-0000057