


Technical Data: Medstrom SoloMH



Height Range	19cm - 81cm
Safe Working Load (SWL)	258kg (41st)
Max Patient Weight (MPW)	193kg (30st)
Back Rest Angle	70°
Knee Angle	30°
Calf Angle	20°
Trend / Reverse Trend Angle	+ 14° / - 14°
Total Length & Width (Side Panels)	225cm x 103cm
Total Length & Width (Side Rails)	227cm x 106cm
Sleep Deck Dimensions (Side Panels)	206cm x 91cm
Sleep Deck Dimensions (Side Rails)	206cm x 93cm
Total Bed Weight (Side Panels)	139kg (22st)
Total Bed Weight (Side Rails)	149kg (23st)



Scan the QR code to access the clinically focused training video